





AUTOMATED ENDOSCOPE REPROCESSOR (AER)
FILTER LIST

OEM MODEL NUMBER	CLINICAL CHOICE® MODEL NUMBER	DESCRIPTION
Compatible with Medivators®		
MF01-0014 	MFO-100	1 micron filter for DSD EDGE & ADVANTAGE type processors
CF01-0017 	MFO-150	1 micron pre-filter for 5000, MV, CLM & CER type processors
MF01-0015 	MFO-200	0.45 micron filter for DSD EDGE & ADVANTAGE type processors
MF-01-003 	MFO-250	5 micron inline pre-sediment water filter for DSD type processors
MF01-0011 	MFO-300	Disinfectant filter (12 per package) for DSD & SSD type processors
MF01-0028 	MFO-375	Air filter (8 per package) for DSD EDGE, ADVANTAGE & ADVANTAGE PLUS type processors
MF01-0013 	MFO-500	0.2 micron bacterial retention grade filter for DSD EDGE & ADVANTAGE type processors
CF01-0018 	MFO-600	0.2 micron post filter for 5000, MV, CLM & CER type processors

OEM MODEL NUMBER	CLINICAL CHOICE® MODEL NUMBER	DESCRIPTION
Compatible with Medivators®		
MF01-0070 	ADP-100	1 micron nominal pre-filter for ADVANTAGE PLUS type processors
MF01-0071 	ADP-200	0.4 micron absolute pre-filter for ADVANTAGE PLUS type processors
MF01-0069 	ADP-500	0.1 micron bacterial retentive pre-filter for ADVANTAGE PLUS type processors
Compatible with Steris®		
A1501 	MF-STE A	« A » pre-filter
A1504 	MF-STE B	« B » pre-filter
A1562 	MF-STE B1E	« B » pre-filter for SYSTEM 1E type processors
A1503 	MF-STE D	Internal air filter

OEM MODEL NUMBER	CLINICAL CHOICE® MODEL NUMBER	DESCRIPTION
Compatible with Olympus®		
MF01-0014PL 	OLY-100	10 inch 1 micron pre-filter for OER Pro®
MF01-0015PL 	OLY-200	10 inch 0.45 micron pre-filter for OER Pro®
MAJ-823 	OLY-300	0.2um Internal Air Filter for OER Pro®
MAJ-824 	OLY-400	Dual Layer 0.65um/0.2um bacterial retention filter with centering ring for OER Pro®
Compatible with J&J® /ASP® AER		
20303 	MFJ-100	1 micron pre-sediment filter (30 per package)
20304 	MFJ-200	0.2 micron bacteria filter
Compatible with Custom Ultrasonics®		
16030 	CU-50	Water and/or air adapter insert (50 per package)
13015 	CU-100	25 micron pre-sediment filter (6 per package)
13016 	CU-200	5 micron pre-sediment filter (6 per package)

OEM MODEL NUMBER	CLINICAL CHOICE® MODEL NUMBER	DESCRIPTION
13017 	CU-300	0.1 micron bacteria filter
21017 	CU-400	20" charcoal replacement filter for SYSTEM 93 type processors
Compatible with Balston®		
200-75-R 	MF-ST200	20" filter cartridges for steam sterilizers (15 per package) *includes 2 rubber gaskets

DISCLAIMER:

Product names, descriptions, and marks that appear in this document may be trademarks and /or registered trademarks of their respective owners, and are used herein solely for informational purposes.